

Product Name: NUCB2 Mouse Monoclonal Antibody**Catalog #: AMM83180**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB,IHC,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	WB 1:500-1:1000,IHC 1:100-1:200,IP 1:100-1:200
Molecular Weight	50.2kDa

Antigen Information


Gene Name	NUCB2
Alternative Names	NUCB2
Gene ID	
SwissProt ID	P80303
Immunogen	Purified recombinant fragment of human NUCB2

Background**Research Area**


Image Data



Western blot analysis using NUCB2 mouse mAb against HL-60 lysate



Lane 1: HL-60 lysate ; Lane 2: NUCB2 immunoprecipitated from HL-60 lysate by AMM83180 ; Lane 3: The same as Lane 2 but use IgG isotype control antibody ;Result: AMM83180 can immunoprecipitate NUCB2;



IHC-P analysis of stomach tissue by NUCB2 antibody (AMM83180)