
Product Name: RGS8 Mouse Monoclonal Antibody**Catalog #: AMM83188**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	IHC
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	IHC 1:100-1:200
Molecular Weight	/

Antigen Information

Gene Name	RGS8
Alternative Names	RGS8
Gene ID	
SwissProt ID	P57771
Immunogen	Purified recombinant fragment of human RGS8

Background**Research Area**

Image Data



IHC-P analysis of cerebellum tissue by RGS8 antibody (RGS8 Mouse Monoclonal Antibody)