

Product Name: CDC25A Mouse Monoclonal Antibody
Catalog #: AMM83189



Summary

Production Name	CDC25A Mouse Monoclonal Antibody
Description	Mouse Monoclonal Antibody
Host	Mouse
Application	IHC
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Immunogen

Gene Name	CDC25A
Alternative Names	CDC25A
Gene ID	
SwissProt ID	P30304. Purified recombinant fragment of human CDC25A

Application

Dilution Ratio	IHC:1:100-1:200
Molecular Weight	/


Background

Product Name: CDC25A Mouse Monoclonal Antibody
Catalog #: AMM83189



Research Area

Image Data



IHC-P analysis of breast cancer tissue by CDC25A antibody (AMM83189)

Note

For research use only.