

Product Name: SERPINC1 Mouse Monoclonal Antibody
Catalog #: AMM83196



Summary

Production Name	SERPINC1 Mouse Monoclonal Antibody
Description	Mouse Monoclonal Antibody
Host	Mouse
Application	IP
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Immunogen

Gene Name	SERPINC1
Alternative Names	SERPINC1
Gene ID	
SwissProt ID	P01008. Purified recombinant fragment of human SERPINC1

Application

Dilution Ratio	IP:1:100-1:200
Molecular Weight	/


Background

Product Name: SERPINC1 Mouse Monoclonal Antibody
Catalog #: AMM83196



Research Area

Image Data



Lane 1: plasma ; Lane 2: SERPINC1 immunoprecipitated from plasma by AMM83196 ; Lane 3: The same as Lane 2 but use IgG isotype control antibody ;Result: AMM83196 can immunoprecipitate SERPINC1;

Note

For research use only.