
Product Name: CCL16 Mouse Monoclonal Antibody**Catalog #: AMM83198**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	WB 1:500-1:1000
Molecular Weight	13.6kDa

Antigen Information

Gene Name	CCL16
Alternative Names	CCL16
Gene ID	
SwissProt ID	O15467
Immunogen	Purified recombinant fragment of human CCL16

Background**Research Area**

Image Data



Result: CCL16 Mouse Monoclonal Antibody can detect human CCL16 by Western blotting.