
Product Name: ABRAXAS2 Mouse Monoclonal Antibody**Catalog #: AMM83204**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	WB 1:500-1:1000,IP 1:100-1:200
Molecular Weight	46.9kDa

Antigen Information


Gene Name	ABRAXAS2
Alternative Names	ABRAXAS2
Gene ID	
SwissProt ID	Q15018
Immunogen	Purified recombinant fragment of human ABRAXAS2

Background**Research Area**

Image Data



Western blot analysis using ABRAXAS2 mouse mAb against Hep G2 lysate



Lane 1: Hep G2 lysate ; Lane 2: ABRAXAS2 immunoprecipitated from Hep G2 lysate by AMM83204 ; Lane 3: The same as Lane 2 but use IgG isotype control antibody ;Result: AMM83204 can immunoprecipitate ABRAXAS2;