
Product Name: WWP1 Mouse Monoclonal Antibody**Catalog #: AMM83206**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	WB 1:500-1:1000
Molecular Weight	105.2kDa

Antigen Information

Gene Name	WWP1
Alternative Names	WWP1
Gene ID	
SwissProt ID	Q9H0M0
Immunogen	Purified recombinant fragment of human WWP1

Background**Research Area**

Image Data



Western blot analysis using WWP1 mouse mAb against HeLa lysate