
Product Name: CMPK2 Mouse Monoclonal Antibody**Catalog #: AMM83211**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB,IHC
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	WB 1:500-1:1000,IHC 1:100-1:200
Molecular Weight	49.4kDa

Antigen Information


Gene Name	CMPK2
Alternative Names	CMPK2
Gene ID	
SwissProt ID	Q5EBM0
Immunogen	Purified recombinant fragment of human CMPK2

Background**Research Area**

Image Data



Western blot analysis using CMPK2 mouse mAb against Jurkat lysate



IHC-P analysis of adrenal gland tissue by CMPK2 antibody (AMM83211)