

Product Name: UBE2Z Mouse Monoclonal Antibody**Catalog #: AMM83213**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB,IHC,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG2b
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	WB 1:500-1:1000,IHC 1:100-1:200,IP 1:100-1:200
Molecular Weight	38.2kDa

Antigen Information


Gene Name	UBE2Z
Alternative Names	UBE2Z
Gene ID	
SwissProt ID	Q9H832
Immunogen	Purified recombinant fragment of human UBE2Z

Background**Research Area**


Image Data



Western blot analysis using UBE2Z mouse mAb against MCF-7 lysate



Lane 1: MCF7 lysate ; Lane 2: UBE2Z immunoprecipitated from MCF7 lysate by AMM83213 ; Lane 3: The same as Lane 2 but use IgG isotype control antibody ;Result: AMM83213 can immunoprecipitate UBE2Z;



IHC-P analysis of testis tissue by UBE2Z antibody (AMM83213)