

Product Name: NUB1 Mouse Monoclonal Antibody**Catalog #: AMM83216**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB,IHC
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	WB 1:500-1:1000,IHC 1:100-1:200
Molecular Weight	70kDa

Antigen Information


Gene Name	NUB1
Alternative Names	NUB1
Gene ID	
SwissProt ID	Q9Y5A7
Immunogen	Purified recombinant fragment of human NUB1

Background**Research Area**

Image Data



Lane 1: 15 μ g of Raji lysate; Lane 2: 15 μ g of SH-SY5Y lysate; Lane 3: 15 μ g of HEK-293 lysate; Result: NUB1 Mouse Monoclonal Antibody can detect NUB1 by Western blotting;



IHC-P analysis of small intestine tissue by NUB1 antibody (AMM83216)