
Product Name: HIP1R Mouse Monoclonal Antibody**Catalog #: AMM83218**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB,IHC,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	WB 1:500-1:1000,IHC 1:100-1:200,IP 1:100-1:200
Molecular Weight	119.4kDa

Antigen Information


Gene Name	HIP1R
Alternative Names	HIP1R
Gene ID	
SwissProt ID	O75146
Immunogen	Purified recombinant fragment of human HIP1R

Background**Research Area**


Image Data



Western blot analysis using HIP1R mouse mAb against A-431 lysate



Lane 1: A-431 lysate ; Lane 2: HIP1R immunoprecipitated from A-431 lysate by AMM83218 ; Lane 3: The same as Lane 2 but use IgG isotype control antibody ;Result: AMM83218 can immunoprecipitate HIP1R;



IHC-P analysis of kidney tissue by HIP1R antibody (AMM83218)