
Product Name: RNF135 Mouse Monoclonal Antibody**Catalog #: AMM83220**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	IHC
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	IHC 1:200-1:500
Molecular Weight	/

Antigen Information

Gene Name	RNF135
Alternative Names	RNF135
Gene ID	
SwissProt ID	Q8IUD6
Immunogen	Purified recombinant fragment of human RNF135

Background**Research Area**

Image Data



IHC-P analysis of kidney tissue by RNF135 antibody (RNF135 Mouse Monoclonal Antibody)