
Product Name: CMTM6 Mouse Monoclonal Antibody**Catalog #: AMM83223**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB,IHC,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	WB 1:500-1:1000,IHC 1:200-1:500,IP 1:100-1:200
Molecular Weight	20kDa

Antigen Information


Gene Name	CMTM6
Alternative Names	CMTM6
Gene ID	
SwissProt ID	Q9NX76
Immunogen	Purified recombinant fragment of human CMTM6

Background**Research Area**


Image Data



Western blot analysis using CMTM6 mouse mAb against Jurkat lysate



Lane 1: Jurkat lysate; Lane 2: CMTM6 immunoprecipitated from Jurkat lysate by AMM83223; Lane 3: The same as Lane 2 but use IgG isotype control antibody; Result: AMM83223 can immunoprecipitate CMTM6;



IHC-P analysis of endometrium tissue by CMTM6 antibody (AMM83223)