
Product Name: MTHFD2 Mouse Monoclonal Antibody**Catalog #: AMM83224**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB,IHC,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	WB 1:500-1:1000,IHC 1:100-1:200,IP 1:100-1:200
Molecular Weight	37.9kDa

Antigen Information


Gene Name	MTHFD2
Alternative Names	MTHFD2
Gene ID	
SwissProt ID	P13995
Immunogen	Purified recombinant fragment of human MTHFD2

Background**Research Area**


Image Data



Western blot analysis using MTHFD2 mouse mAb against A-431 lysate



Lane 1: A-431 lysate ; Lane 2: MTHFD2 immunoprecipitated from A-431 lysate by AMM83224 ; Lane 3: The same as Lane 2 but use IgG isotype control antibody ;Result: AMM83224 can immunoprecipitate MTHFD2;



IHC-P analysis of kidney tissue by MTHFD2 antibody (AMM83224)