
Product Name: GSTT2B Mouse Monoclonal Antibody**Catalog #: AMM83225**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB
Reactivity	Human, Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG2a
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000
Molecular Weight	27.5kDa

Antigen Information

Gene Name	GSTT2B
Alternative Names	GSTT2B
Gene ID	
SwissProt ID	P0CG30
Immunogen	Purified recombinant fragment of human GSTT2B

Background**Research Area**

Image Data



Western blot analysis using GSTT2B mouse mAb against rat liver tissue lysate