
Product Name: ZFP36 Mouse Monoclonal Antibody**Catalog #: AMM83227**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG2a
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	WB 1:500-1:1000,IP 1:100-1:200
Molecular Weight	34kDa

Antigen Information


Gene Name	ZFP36
Alternative Names	ZFP36
Gene ID	
SwissProt ID	P26651
Immunogen	Purified recombinant fragment of human ZFP36

Background**Research Area**

Image Data



Western blot analysis using ZFP36 mouse mAb against HCT116 lysate



Lane 1: HCT116 lysate ; Lane 2: ZFP36 immunoprecipitated from HCT116 lysate by AMM83227 ; Lane 3: The same as Lane 2 but use IgG isotype control antibody ;Result: AMM83227 can immunoprecipitate ZFP36;