

Product Name: CXCL12 Mouse Monoclonal Antibody**Catalog #: AMM83232**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB,IHC
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG2a
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	WB 1:500-1:1000,IHC 1:100-1:200
Molecular Weight	11kDa

Antigen Information


Gene Name	CXCL12
Alternative Names	CXCL12
Gene ID	
SwissProt ID	P48061
Immunogen	Purified recombinant fragment of human CXCL12

Background**Research Area**

Image Data



Result: CXCL12 Mouse Monoclonal Antibody can detect CXCL12 by Western blotting.



IHC-P analysis of tonsil tissue by CXCL12 antibody (AMM83232)