

Product Name: Mitf (phospho S157) Mouse Monoclonal Antibody**Catalog #: AMM83238**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	IHC
Reactivity	Rat
Conjugation	Unconjugated
Modification	Phosphorylated
Isotype	Mouse IgG2b
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	IHC 1:50-1:200
Molecular Weight	/

Antigen Information

Gene Name	Mitf (phospho S157)
Alternative Names	Mitf (phospho S157)
Gene ID	
SwissProt ID	O88368
Immunogen	Mitf peptide

Background**Research Area**

Image Data



IHC-P analysis of cutaneous melanoma tissue by Mitf antibody (Mitf (phospho S157) Mouse Monoclonal Antibody)