

Product Name: UBR5 Mouse Monoclonal Antibody
Catalog #: AMM83240



Summary

Production Name	UBR5 Mouse Monoclonal Antibody
Description	Mouse Monoclonal Antibody
Host	Mouse
Application	IHC
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG2a
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Immunogen

Gene Name	UBR5
Alternative Names	UBR5
Gene ID	
SwissProt ID	O95071. Purified recombinant fragment of human UBR5

Application

Dilution Ratio	IHC:1:100-1:200
Molecular Weight	/


Background

Product Name: UBR5 Mouse Monoclonal Antibody
Catalog #: AMM83240



Research Area

Image Data



IHC-P analysis of testis tissue by UBR5 antibody (AMM83240)

Note

For research use only.