

Product Name: TNFRSF10C Mouse Monoclonal Antibody
Catalog #: AMM83241

Summary

Production Name	TNFRSF10C Mouse Monoclonal Antibody
Description	Mouse Monoclonal Antibody
Host	Mouse
Application	IHC
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG2b
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Immunogen

Gene Name	TNFRSF10C
Alternative Names	TNFRSF10C
Gene ID	
SwissProt ID	O14798. Purified recombinant fragment of human TNFRSF10C

Application

Dilution Ratio	IHC:1:500-1:1000
Molecular Weight	/


Background

Product Name: TNFRSF10C Mouse Monoclonal Antibody
Catalog #: AMM83241



Research Area

Image Data



IHC-P analysis of spleen tissue by TNFRSF10C antibody (AMM83241)

Note

For research use only.