
Product Name: CAPS Mouse Monoclonal Antibody**Catalog #: AMM83245**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB,IHC,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG2a
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	WB 1:500-1:1000,IHC 1:500-1:1000,IP 1:100-1:200
Molecular Weight	30.2kDa

Antigen Information


Gene Name	CAPS
Alternative Names	CAPS
Gene ID	
SwissProt ID	Q13938
Immunogen	Purified recombinant fragment of human CAPS

Background**Research Area**


Image Data



Western blot analysis using CAPS mouse mAb against U-251MG lysate



Lane 1: U-251MG lysate ; Lane 2: CAPS immunoprecipitated from U-251MG lysate by AMM83245 ; Lane 3: The same as Lane 2 but use IgG isotype control antibody ;Result: AMM83245 can immunoprecipitate CAPS;



IHC-P analysis of fallopian tube tissue by CAPS antibody (AMM83245)