

Product Name: SH3RF1 Mouse Monoclonal Antibody**Catalog #: AMM83246**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB,IP
Reactivity	Human, Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG2a
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	WB 1:500-1:1000,IP 1:100-1:200
Molecular Weight	93.1kDa

Antigen Information


Gene Name	SH3RF1
Alternative Names	SH3RF1
Gene ID	
SwissProt ID	Q7Z6J0
Immunogen	Purified recombinant fragment of human SH3RF1

Background**Research Area**

Image Data



Western blot analysis using SH3RF1 mouse mAb against HeLa lysate and mouse brain tissue lysate



Lane 1: Mouse brain tissue lysate ; Lane 2: SH3RF1 immunoprecipitated from mouse brain tissue lysate by AMM83246 ; Lane 3: The same as Lane 2 but use IgG isotype control antibody ;Result: AMM83246 can immunoprecipitate SH3RF1;