

Product Name: MAP3K1 Mouse Monoclonal Antibody**Catalog #: AMM83250**

For research use only.

Summary

| | |
|----------------------|---|
| Description | Mouse monoclonal Antibody |
| Host | Mouse |
| Application | IHC |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Modification | Unmodified |
| Isotype | Mouse IgG1 |
| Clonality | Monoclonal |
| Form | Liquid |
| Concentration | 1mg/ml |
| Storage | Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles. |
| Shipping | Ice bags |
| Buffer | Purified antibody in PBS with 0.05% sodium azide |
| Purification | Affinity Purification |

Application


| | |
|-------------------------|-----------------|
| Dilution Ratio | IHC 1:100-1:200 |
| Molecular Weight | / |

Antigen Information

| | |
|--------------------------|---|
| Gene Name | MAP3K1 |
| Alternative Names | MAP3K1 |
| Gene ID | |
| SwissProt ID | Q13233 |
| Immunogen | Purified recombinant fragment of human MAP3K1 |

Background**Research Area**

Image Data



IHC-P analysis of liver tissue by MAP3K1 antibody (MAP3K1 Mouse Monoclonal Antibody)