
Product Name: TRIML2 Mouse Monoclonal Antibody**Catalog #: AMM83254**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	IHC,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG2a
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	IHC 1:100-1:200,IP 1:100-1:200
Molecular Weight	/

Antigen Information


Gene Name	TRIML2
Alternative Names	TRIML2
Gene ID	
SwissProt ID	Q8N7C3
Immunogen	Purified recombinant fragment of human TRIML2

Background**Research Area**

Image Data



Lane 1: A549 lysate ; Lane 2: TRIML2 immunoprecipitated from A549 lysate by K52045_1P15007 ; Lane 3: The same Lane 2 but use IgG isotype control antibody ;Result: K52045_1P15007 can immunoprecipitate TRIML2;



IHC-P analysis of testis tissue by TRIML2 antibody (K52045_1P15007)