
Product Name: Acetyl-SMC2 (Lys114) Mouse Monoclonal Antibody**Catalog #: AMM83264**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Acetylated
Isotype	Mouse IgG2a
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:1000,IP 1:100-1:200
Molecular Weight	135.7kDa

Antigen Information


Gene Name	SMC2 (acetyl K114)
Alternative Names	SMC2 (acetyl K114)
Gene ID	
SwissProt ID	O95347
Immunogen	Human SMC2 peptide

Background**Research Area**

Image Data



Western blot analysis using SMC2 mouse mAb against Raji lysate



Lane 1: Raji lysate; Lane 2: SMC2 immunoprecipitated from Raji lysate by P151963; Lane 3: The same as Lane 2 but use IgG isotype control antibody; Result: P151963 can immunoprecipitate SMC2;