
Product Name: BCL2L15 Mouse Monoclonal Antibody**Catalog #: AMM83266**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB,IHC,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	WB 1:500-1:1000,IHC 1:100-1:200,IP 1:100-1:200
Molecular Weight	17.7kDa

Antigen Information


Gene Name	BCL2L15
Alternative Names	BCL2L15
Gene ID	
SwissProt ID	Q5TBC7
Immunogen	Purified recombinant fragment of human BCL2L15

Background**Research Area**


Image Data



Western blot analysis using BCL2L15 mouse mAb against HT-29 lysate



Lane 1: HT-29 lysate ; Lane 2: BCL2L15 immunoprecipitated from HT-29 lysate by AMM83266 ; Lane 3: The same as Lane 2 but use IgG isotype control antibody ;Result: AMM83266 can immunoprecipitate BCL2L15;



IHC-P analysis of colon tissue by BCL2L15 antibody (AMM83266)