

Product Name: SULT2B1 Mouse Monoclonal Antibody**Catalog #: AMM83272**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	IP 1:100-1:200
Molecular Weight	/

Antigen Information

Gene Name	SULT2B1
Alternative Names	SULT2B1
Gene ID	
SwissProt ID	O00204
Immunogen	Purified recombinant fragment of human SULT2B1

Background**Research Area**

Image Data



Lane 1: HT-29 lysate ; Lane 2: SULT2B1 immunoprecipitated from HT-29 lysate by SULT2B1 Mouse Monoclonal Antibody ; Lane 3: The same as Lane 2 but use IgG isotype control antibody ;Result: SULT2B1 Mouse Monoclonal Antibody can immunoprecipitate SULT2B1;