
Product Name: DTX3L Mouse Monoclonal Antibody**Catalog #: AMM83276**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	WB 1:500-1:1000,IP 1:100-1:200
Molecular Weight	83.6kDa

Antigen Information


Gene Name	DTX3L
Alternative Names	DTX3L
Gene ID	
SwissProt ID	Q8TDB6
Immunogen	Purified recombinant fragment of human DTX3L

Background**Research Area**

Image Data



Western blot analysis using DTX3L mouse mAb against A-431 lysate



Lane 1: A-431 lysate ; Lane 2: DTX3L immunoprecipitated from A-431 lysate by AMM83276 ; Lane 3: The same as Lane 2 but use IgG isotype control antibody ;Result: AMM83276 can immunoprecipitate DTX3L;