
Product Name: C1QTNF1 Mouse Monoclonal Antibody**Catalog #: AMM83283**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	IHC
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG2b
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	IHC 1:100-1:200
Molecular Weight	/

Antigen Information

Gene Name	C1QTNF1
Alternative Names	C1QTNF1
Gene ID	
SwissProt ID	Q9BXJ1
Immunogen	Purified recombinant fragment of human C1QTNF1

Background**Research Area**

Image Data



IHC-P analysis of placenta tissue by C1QTNF1 antibody (C1QTNF1 Mouse Monoclonal Antibody)