

Product Name: PRKCD Mouse Monoclonal Antibody
Catalog #: AMM83285



Summary

Production Name	PRKCD Mouse Monoclonal Antibody
Description	Mouse Monoclonal Antibody
Host	Mouse
Application	WB,IP,IHC
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Immunogen

Gene Name	PRKCD
Alternative Names	PRKCD
Gene ID	
SwissProt ID	Q05655.Purified recombinant fragment of human PRKCD

Application

Dilution Ratio	WB:1:500-1:1000,IP:1:100-1:200,IHC:1:100-1:200
Molecular Weight	77kDa


Background

Product Name: PRKCD Mouse Monoclonal Antibody
Catalog #: AMM83285




Research Area


Image Data



Western blot analysis using PRKCD mouse mAb against Jurkat lysate



Lane 1: Jurkat lysate ; Lane 2: PRKCD immunoprecipitated from Jurkat lysate by AMM83285 ; Lane 3: The same as Lane 2 but use IgG isotype control antibody ;Result: AMM83285 can immunoprecipitate PRKCD;



IHC-P analysis of fallopian tube tissue by PRKCD antibody (AMM83285)

Product Name: PRKCD Mouse Monoclonal Antibody
Catalog #: AMM83285



Note

For research use only.