

Product Name: IFRD1 Mouse Monoclonal Antibody**Catalog #: AMM83286**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB
Reactivity	Human, Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG2a
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	WB 1:500-1:1000
Molecular Weight	Calculated MW: 50 kDa; Observed MW: 50 kDa

Antigen Information

Gene Name	IFRD1
Alternative Names	IFRD1
Gene ID	
SwissProt ID	O00458
Immunogen	Purified recombinant fragment of human IFRD1

Background**Research Area**

Image Data



Western blot analysis using IFRD1 mouse mAb against mouse skeletal muscle tissue lysate