

Product Name: FH Mouse Monoclonal Antibody**Catalog #: AMM83287**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG2a
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	WB 1:500-1:1000,IP 1:100-1:200
Molecular Weight	50kDa

Antigen Information


Gene Name	FH
Alternative Names	FH
Gene ID	
SwissProt ID	P07954
Immunogen	Purified recombinant fragment of human FH

Background**Research Area**

Image Data



Western blot analysis using FH mouse mAb against Ramos lysate



Lane 1: Ramos lysate ; Lane 2: FH immunoprecipitated from Ramos lysate by AMM83287 ; Lane 3: The same as Lane 2 but use IgG isotype control antibody ;Result: AMM83287 can immunoprecipitate FH;