
Product Name: MMAB Mouse Monoclonal Antibody**Catalog #: AMM83288**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	IHC,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG2b
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	IHC 1:100-1:200,IP 1:100-1:200
Molecular Weight	/

Antigen Information


Gene Name	MMAB
Alternative Names	MMAB
Gene ID	
SwissProt ID	Q96EY8
Immunogen	Purified recombinant fragment of human MMAB

Background**Research Area**

Image Data



Lane 1: A549 lysate ; Lane 2: MMAB immunoprecipitated from A549 lysate by MMAB Mouse Monoclonal Antibody ; Lane 3: The same Lane 2 but use IgG isotype control antibody ;Result: MMAB Mouse Monoclonal Antibody can immunoprecipitate MMAB;



IHC-P analysis of liver tissue by MMAB antibody (AMM83288)