

Product Name: CBLL2 Mouse Monoclonal Antibody
Catalog #: AMM83291



Summary

Production Name	CBLL2 Mouse Monoclonal Antibody
Description	Mouse Monoclonal Antibody
Host	Mouse
Application	WB,IP
Reactivity	Human, Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG2a
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Immunogen

Gene Name	CBLL2
Alternative Names	CBLL2
Gene ID	
SwissProt ID	Q8N7E2.Purified recombinant fragment of human CBLL2

Application

Dilution Ratio	WB:1:500-1:1000,IP:1:100-1:200
Molecular Weight	48.8kDa


Background

Product Name: CBLL2 Mouse Monoclonal Antibody
Catalog #: AMM83291




Research Area

Image Data



Western blot analysis using CBLL2 mouse mAb against mouse testis tissue lysate



Lane 1: Mouse testis tissue lysate ; Lane 2: CBLL2 immunoprecipitated from mouse testis tissue lysate by AMM83291 ; Lane 3: The same as Lane 2 but use IgG isotype control antibody ;Result: AMM83291 can immunoprecipitate CBLL2;

Note

For research use only.