
Product Name: NPTX1 Mouse Monoclonal Antibody**Catalog #: AMM83295**

For research use only.

Summary

| | |
|----------------------|---|
| Description | Mouse monoclonal Antibody |
| Host | Mouse |
| Application | WB,IP |
| Reactivity | Human, Rat |
| Conjugation | Unconjugated |
| Modification | Unmodified |
| Isotype | Mouse IgG1 |
| Clonality | Monoclonal |
| Form | Liquid |
| Concentration | 1mg/ml |
| Storage | Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles. |
| Shipping | Ice bags |
| Buffer | Purified antibody in PBS with 0.05% sodium azide |
| Purification | Affinity Purification |

Application


| | |
|-------------------------|--------------------------------|
| Dilution Ratio | WB 1:500-1:1000,IP 1:100-1:200 |
| Molecular Weight | 47.1kDa |

Antigen Information


| | |
|--------------------------|--|
| Gene Name | NPTX1 |
| Alternative Names | NPTX1 |
| Gene ID | |
| SwissProt ID | Q15818 |
| Immunogen | Purified recombinant fragment of human NPTX1 |

Background**Research Area**

Image Data



Western blot analysis using NPTX1 mouse mAb against rat brain



Lane 1: Rat brain tissue lysate ; Lane 2: NPTX1 immunoprecipitated from rat brain tissue lysate by AMM83295 ; Lane 3: The same as Lane 2 but use IgG isotype control antibody ;Result: AMM83295 can immunoprecipitate NPTX1;