
Product Name: FGF20 Mouse Monoclonal Antibody**Catalog #: AMM83296**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG2b
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	WB 1:500-1:1000
Molecular Weight	Calculated MW: 23 kDa; Observed MW: 23 kDa

Antigen Information

Gene Name	FGF20
Alternative Names	FGF20
Gene ID	
SwissProt ID	Q9NP95
Immunogen	Purified recombinant fragment of human FGF20

Background**Research Area**

Image Data



Result: FGF20 Mouse Monoclonal Antibody can detect FGF20 by Western blotting.