

Product Name: ZNF598 Mouse Monoclonal Antibody
Catalog #: AMM83297



Summary

Production Name	ZNF598 Mouse Monoclonal Antibody
Description	Mouse Monoclonal Antibody
Host	Mouse
Application	WB,IP
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Immunogen

Gene Name	ZNF598
Alternative Names	ZNF598
Gene ID	
SwissProt ID	Q86UK7.Purified recombinant fragment of human ZNF598

Application

Dilution Ratio	WB:1:500-1:1000,IP:1:100-1:200
Molecular Weight	99kDa


Background

Product Name: ZNF598 Mouse Monoclonal Antibody
Catalog #: AMM83297




Research Area

Image Data



Western blot analysis using ZNF598 mouse mAb against Jurkat lysate



Lane 1: Jurkat lysate ; Lane 2: ZNF598 immunoprecipitated from Jurkat lysate by AMM83297 ; Lane 3: The same as Lane 2 but use IgG isotype control antibody ;Result: AMM83297 can immunoprecipitate ZNF598;

Note

For research use only.