

Product Name: TOMM20 Mouse Monoclonal Antibody
Catalog #: AMM83301



Summary

Production Name	TOMM20 Mouse Monoclonal Antibody
Description	Mouse Monoclonal Antibody
Host	Mouse
Application	WB,IP,IHC
Reactivity	Human, Rat, Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Immunogen

Gene Name	TOMM20
Alternative Names	TOMM20
Gene ID	
SwissProt ID	Q15388.Purified recombinant fragment of human TOMM20

Application

Dilution Ratio	WB:1:500-1:1000,IP:1:100-1:200,IHC:1:100-1:200
Molecular Weight	Calculated MW: 16 kDa; Observed MW: 16 kDa


Background

Product Name: TOMM20 Mouse Monoclonal Antibody
Catalog #: AMM83301




Research Area


Image Data



Western blot analysis using TOMM20 mouse mAb against mouse brain and rat brain



Western blot analysis using TOMM20 mouse mAb against HEK-293 lysate



IHC-P analysis of kidney tissue by TOMM20 antibody (AMM83301)

Note

Product Name: TOMM20 Mouse Monoclonal Antibody
Catalog #: AMM83301



For research use only.