

Product Name: CCL8 Mouse Monoclonal Antibody
Catalog #: AMM83307



Summary

Production Name	CCL8 Mouse Monoclonal Antibody
Description	Mouse Monoclonal Antibody
Host	Mouse
Application	IHC
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Immunogen

Gene Name	CCL8
Alternative Names	CCL8
Gene ID	
SwissProt ID	P80075. Purified recombinant fragment of human CCL8

Application

Dilution Ratio	IHC:1:100-1:200
Molecular Weight	/


Background

Product Name: CCL8 Mouse Monoclonal Antibody
Catalog #: AMM83307



Research Area

Image Data



IHC-P analysis of small intestine tissue by CCL8 antibody (AMM83307)

Note

For research use only.