

**Product Name: SARS-CoV-2 Nucleocapsid Protein
Mouse Monoclonal Antibody
Catalog #: AMM83308**

Summary

Production Name	SARS-CoV-2 Nucleocapsid Protein Mouse Monoclonal Antibody
Description	Mouse Monoclonal Antibody
Host	Mouse
Application	WB
Reactivity	SARS-CoV-2

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG2b
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Immunogen

Gene Name	SARS-CoV-2 Nucleocapsid Protein
Alternative Names	SARS-CoV-2 Nucleocapsid Protein
Gene ID	
SwissProt ID	P0DTC9. Recombinant SARS-CoV-2 nucleocapsid protein

Application


Dilution Ratio	WB:1:500-1:1000
Molecular Weight	55kDa

Background

Product Name: SARS-CoV-2 Nucleocapsid Protein
Mouse Monoclonal Antibody
Catalog #: AMM83308

Research Area

Image Data



Result: AMM83308 can detect recombinant nucleocapsid protein by Western blotting.

Note

For research use only.