
Product Name: RNF146 Mouse Monoclonal Antibody**Catalog #: AMM83312**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	IP
Reactivity	Human, Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG2b
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	IP 1:100-1:200
Molecular Weight	/

Antigen Information

Gene Name	RNF146
Alternative Names	RNF146
Gene ID	
SwissProt ID	Q9NTX7
Immunogen	Purified recombinant fragment of human RNF146

Background**Research Area**

Image Data



Lane 1: Rat brain lysate ; Lane 2: RNF146 immunoprecipitated from rat brain lysate by RNF146 Mouse Monoclonal Antibody ; Lane 3: The same as Lane 2 but use IgG isotype control antibody ;Result: RNF146 Mouse Monoclonal Antibody can immunoprecipitate RNF146;