

Product Name: DPCD Mouse Monoclonal Antibody**Catalog #: AMM83315**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB,IP
Reactivity	Human, Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG2b
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	WB 1:500-1:1000,IP 1:100-1:200
Molecular Weight	23.2kDa

Antigen Information


Gene Name	DPCD
Alternative Names	DPCD
Gene ID	
SwissProt ID	Q9BVM2
Immunogen	Purified recombinant fragment of human DPCD

Background**Research Area**

Image Data



Western blot analysis using DPCD mouse mAb against rat testis lysate



Lane 1: Mouse testis tissue lysate ; Lane 2: DPCD immunoprecipitated from mouse testis tissue lysate by AMM83315 ; Lane 3: The same as Lane 2 but use IgG isotype control antibody ;Result: AMM83315 can immunoprecipitate DPCD;