
Product Name: WDR59 Mouse Monoclonal Antibody**Catalog #: AMM83316**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG2a
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	WB 1:1000-1:2000
Molecular Weight	109.8kDa

Antigen Information

Gene Name	WDR59
Alternative Names	WDR59
Gene ID	
SwissProt ID	Q6PJ19
Immunogen	Purified recombinant fragment of human WDR59

Background**Research Area**

Image Data



Western blot analysis using WDR59 mouse mAb against Jurkat lysate