

Product Name: LRRN1 Mouse Monoclonal Antibody**Catalog #: AMM83317**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	WB 1:500-1:1000,IP 1:100-1:200
Molecular Weight	80.7kDa

Antigen Information


Gene Name	LRRN1
Alternative Names	LRRN1
Gene ID	
SwissProt ID	Q6UXK5
Immunogen	Purified recombinant fragment of human LRRN1

Background**Research Area**

Image Data



Western blot analysis using LRRN1 mouse mAb against SH-SY5Y lysate and brain lysate



Lane 1: SH-SY5Y lysate ; Lane 2: LRRN1 immunoprecipitated from SH-SY5Y lysate by AMM83317 ; Lane 3: The same as Lane 2 but use IgG isotype control antibody ;Result: AMM83317 can immunoprecipitate LRRN1;