

Product Name: COLEC12 Mouse Monoclonal Antibody
Catalog #: AMM83320



Summary

Production Name	COLEC12 Mouse Monoclonal Antibody
Description	Mouse Monoclonal Antibody
Host	Mouse
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG2a
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Immunogen

Gene Name	COLEC12
Alternative Names	COLEC12
Gene ID	
SwissProt ID	Q5KU26. Purified recombinant fragment of human COLEC12

Application

Dilution Ratio	WB:1:500-1:1000
Molecular Weight	81.5kDa


Background

Product Name: COLEC12 Mouse Monoclonal Antibody
Catalog #: AMM83320



Research Area

Image Data



Western blot analysis using COLEC12 mouse mAb against Raji lysate

Note

For research use only.