
Product Name: NAGK Mouse Monoclonal Antibody**Catalog #: AMM83325**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB,IHC,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG2a
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	WB 1:500-1:1000,IHC 1:100-1:200,IP 1:100-1:200
Molecular Weight	37.4kDa

Antigen Information


Gene Name	NAGK
Alternative Names	NAGK
Gene ID	
SwissProt ID	Q9UJ70
Immunogen	Purified recombinant fragment of human NAGK

Background**Research Area**


Image Data



Lane 1: 15 μ g of HeLa lysate; Lane 2: 15 μ g of Jurkat lysate; Lane 3: 15 μ g of HEK-293 lysate; Lane 4: 15 μ g of A-431 lysate; Result: NAGK Mouse Monoclonal Antibody can detect NAGK by Western blotting.



Lane 1: HeLa lysate; Lane 2: NAGK immunoprecipitated from HeLa lysate by AMM83325; Lane 3: The same as Lane 2 but use IgG isotype control antibody; Result: AMM83325 can immunoprecipitate NAGK;



IHC-P analysis of tonsil tissue by NAGK antibody (AMM83325)