
Product Name: CHCHD6 Mouse Monoclonal Antibody**Catalog #: AMM83329**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	IHC,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG2a
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	IHC 1:100-1:200,IP 1:100-1:200
Molecular Weight	/

Antigen Information


Gene Name	CHCHD6
Alternative Names	CHCHD6
Gene ID	
SwissProt ID	Q9BRQ6
Immunogen	Purified recombinant fragment of human CHCHD6

Background**Research Area**

Image Data



Lane 1: Jurkat lysate ; Lane 2: CHCHD6 immunoprecipitated from Jurkat lysate by CHCHD6 Mouse Monoclonal Antibody ; Lane 3: The same Lane 2 but use IgG isotype control antibody ;Result: CHCHD6 Mouse Monoclonal Antibody can immunoprecipitate CHCHD6;



IHC-P analysis of stomach tissue by CHCHD6 antibody (AMM83329)