

**Product Name: DECR1 Mouse Monoclonal Antibody****Catalog #: AMM83332**

For research use only.

**Summary**

|                      |                                                                             |
|----------------------|-----------------------------------------------------------------------------|
| <b>Description</b>   | Mouse monoclonal Antibody                                                   |
| <b>Host</b>          | Mouse                                                                       |
| <b>Application</b>   | IHC,IP                                                                      |
| <b>Reactivity</b>    | Human                                                                       |
| <b>Conjugation</b>   | Unconjugated                                                                |
| <b>Modification</b>  | Unmodified                                                                  |
| <b>Isotype</b>       | Mouse IgG1                                                                  |
| <b>Clonality</b>     | Monoclonal                                                                  |
| <b>Form</b>          | Liquid                                                                      |
| <b>Concentration</b> | 1mg/ml                                                                      |
| <b>Storage</b>       | Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles. |
| <b>Shipping</b>      | Ice bags                                                                    |
| <b>Buffer</b>        | Purified antibody in PBS with 0.05% sodium azide                            |
| <b>Purification</b>  | Affinity Purification                                                       |

**Application**


|                         |                                |
|-------------------------|--------------------------------|
| <b>Dilution Ratio</b>   | IHC 1:100-1:200,IP 1:100-1:200 |
| <b>Molecular Weight</b> | /                              |

**Antigen Information**


|                          |                                              |
|--------------------------|----------------------------------------------|
| <b>Gene Name</b>         | DECR1                                        |
| <b>Alternative Names</b> | DECR1                                        |
| <b>Gene ID</b>           |                                              |
| <b>SwissProt ID</b>      | Q16698                                       |
| <b>Immunogen</b>         | Purified recombinant fragment of human DECR1 |

**Background****Research Area**

## Image Data



Lane 1: SK-BR-3 lysate; Lane 2: DECR1 immunoprecipitated from SK-BR-3 by DECR1 Mouse Monoclonal Antibody; Lane 3: The same Lane 2 but use IgG isotype control antibody ;Result: DECR1 Mouse Monoclonal Antibody can immunoprecipitate DECR1;



IHC-P analysis of liver tissue by DECR1 antibody (AMM83332)