
Product Name: IRF9 Mouse Monoclonal Antibody**Catalog #: AMM83337**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	IHC,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	IHC 1:100-1:200,IP 1:100-1:200
Molecular Weight	/

Antigen Information


Gene Name	IRF9
Alternative Names	IRF9
Gene ID	
SwissProt ID	Q00978
Immunogen	Purified recombinant fragment of human IRF9

Background**Research Area**

Image Data



Lane 1: THP-1 lysate ; Lane 2: IRF9 immunoprecipitated from THP-1 lysate by IRF9 Mouse Monoclonal Antibody ; Lane 3: The same Lane 2 but use IgG isotype control antibody ;Result: IRF9 Mouse Monoclonal Antibody can immunoprecipitate IRF9;



IHC-P analysis of lung tissue by IRF9 antibody (AMM83337)