
Product Name: ARFGAP3 Mouse Monoclonal Antibody**Catalog #: AMM83342**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	IHC,IP
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG2b
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Application


Dilution Ratio	IHC 1:500-1:1000,IP 1:100-1:200
Molecular Weight	/

Antigen Information


Gene Name	ARFGAP3
Alternative Names	ARFGAP3
Gene ID	
SwissProt ID	Q9NP61
Immunogen	Purified recombinant fragment of human ARFGAP3

Background**Research Area**

Image Data



Lane 1: Ramos lysate; Lane 2: ARFGAP3 immunoprecipitated from Ramos by ARFGAP3 Mouse Monoclonal Antibody; Lane 3: The same Lane 2 but use IgG isotype control antibody ;Result: ARFGAP3 Mouse Monoclonal Antibody can immunoprecipitate ARFGAP3;



IHC-P analysis of salivary gland tissue by ARFGAP3 antibody (AMM83342)